

NOVASONIX TECHNOLOGY

Tanit is an OUTSIDE IN body contouring platform that combines dual laser stimulation with the CBA current (bioactive current) and the exclusive Slaag Motion 400 system.

OUTSIDE IN technology is a unique system developed by Novasonix, aimed at naturally stimulating and activating the internal cellular metabolism. It is combined with the Slaag Motion 400 rotational massage system, which adds an external mechanical stimulation, thus increasing mobilisation and drainage.

The combination of the three technologies enhances the efficacy of the treatment, achieving visible results from the first session.



Dr. Pablo Naranjo

The Tanit platform has made a difference at my clinic. The combination of laser technologies and mechanical massage achieves truly extraordinary results. My patients are very happy with the changes to their figure.

TANIT PLATFORM

The Tanit machine is designed for medical-aesthetic doctors and specialists wanting to equip their clinic with next-generation Spanish technology.

The results provided by the machine are visible in just a few sessions, creating happy, loyal patients and a high ROI (return on investment) in a short space of time.



NOVASONIX, YOUR TRUSTED PARTNER

Novasonix is a Spanish company with extensive experience in the health and aesthetic medicine sector. With more than 10 years' experience in manufacturing medical-aesthetic devices, we have a presence in over 12 countries, basing our growth on trust, results and professionalism.

Our main aim is to offer maximum innovation in technology, resulting from constant investment in R+D+i, and developing our own patents, such as the CBA current or the NSX method, which are guaranteed by the highest ISO quality certifications.

Our clients are another of the main pillars of our company, which is why we have developed the Novasonix Partners plan. A world of services offered to our clients, providing them with everything from marketing services to ongoing training.



Made
in Spain



High
return



Technical
support



R+D: Own
patented technology



Guaranteed
results



A NEW BEGINNING FOR YOUR SHAPE

*Innovative laser technology that
challenges body contouring*

Tanit

by novasonix

HOW DOES IT WORK?

In its commitment to technological innovation, Novasonix has developed the NSX OUTSIDE IN method. This is a union of technologies aimed at body contouring, drainage, lipolysis, cell regeneration and firming.

NSX METHOD – OUTSIDE TECHNOLOGY

The NSX OUTSIDE method is based on Slaag Motion 400 technology. This is a biomechanical rotation action with a 45-smart sphere head with integrated control. It allows a speed of up to 400 revolutions per minute.

The combination with the Natural Sphere bioactives and their four specific formulations increases the benefits of the treatments performed.



NSX METHOD – IN TECHNOLOGY

The IN method involves dual soft laser emission combined with the CBA current (bioactive current). The double laser wavelength actively promotes the removal of accumulated fats in a natural, non-invasive way. Tanit comprises the following technologies:

- 639nm red light laser emission: stimulates photolipolysis and increases cell metabolism.
- 830nm infrared light: produces adipocyte membrane photo-permeabilization.
- Damped CBA current: induces membrane electroporation.



Tanit

by novasonix

TANIT BENEFITS

Cell biomodulation and increased membrane permeability: photo-stimulation activates the cell mitochondria, increasing cell ATP.

Tissue electroporation: the effect of the current modifies the polarisation of lymphatic channels and vessels, increasing permeability and reducing oedema.

Dermal regeneration and firming: the combination of the dual soft laser and polarised wave stimulates collagen generation and dermal oxygenation.

Adipose and fibrous tissue reduction: the combination of the soft laser and the mechanical effect of the massage acts directly on adipose tissue, reducing and removing fatty acids.

Improved circulation and muscle tone: the 45 smart spheres with integrated control impact on the muscle system, improving its tone, elasticity and firmness.

Cellulite treatment: the union of technologies enhances drainage and microcirculation, as well as deeply mobilising fibrous tissue. This action is a determining factor in the correct treatment of cellulite.

Loss of volume: the synergy of the OUTSIDE IN treatment promotes contouring by reducing the body volume.

NATURAL SPHERE + ACTIVE BODY

The combination of technology and active ingredients acts as a carrier, visibly improving the treatment results. **Natural Sphere** is a range of 4 100% plant-based products with a high active ingredient content. **Active Body** is a treatment with an exclusive concentrated formula combining the most powerful encapsulated slimming active ingredients with an immediate effect.



Tanit

by novasonix

RESULTS

At Novasonix we have over 10 years' experience in engineering and R+D+i, which enables us to develop our own technologies for the medical-aesthetic sector. Our studies have been developed in collaboration with prestigious doctors and health specialists to achieve the best protocols and treatments.



Treated area: Saddlebags
Result: 2cm reduction.
No. sessions: 3



Treated area: Stomach
Result: 3cm reduction.
No. sessions: 1

TECHNICAL FEATURES

TECNOLOGÍA
OUTSIDE

Emission type:	Optical power, 6 pads with 13 units; 10 x 25mw units at 639nm and 3 x 830mw units at 50mw (+0.36mw).
High-voltage pulse current (HVPC) system:	SVPP post-insulation application voltage/O.IA output current.

TECNOLOGÍA
IN

Input voltage / frequency:	100-240v. / 50-60Hz.
Spheres / Maximum speed:	45 Teflon spheres. / 275 rpm.